


Date:	09.12.2022								
Editor:	Füllmann, Robert								
Material:	Scale plate (oxide) from heat treatment furnace from base material 18CrNiMo6-7								
Mounting:	Cold, D=40 mm, EpoClear (epoxy resin) with post-infiltration								
<b>TIP // For reference images, take a look at our wide-ranging microstructure database. There you will find many micrographs of various materials, including their state of heat treatment and the etching agent we used.</b>									
<a href="#">Link to our Microstructure Database (Gefügedatenbank)</a>									
	Grinding					Polishing			
	1	2	3	4	5	1	2	Kurth	
Base	SiC paper	SiC paper	SiC paper	SiC paper	SiC paper	PanTec	PanTec	MolTec	
Grit size	P120	P180	P800	P2500	P4000	6 µm mkd water-based	3 µm mkd water-based	1 µm mkd water-based	
Lubricant	H <sub>2</sub> O	H <sub>2</sub> O	H <sub>2</sub> O	H <sub>2</sub> O	H <sub>2</sub> O	coolTec I	coolTec I	coolTec I	
Drive	150 rpm	150 rpm	150 rpm	150 rpm	150 rpm	150 rpm	150 rpm	50 rpm	
Rotation direction	>>	>>	>>	>>	>>	<<	<<	oscillating	
Downforce !centralpressure!	30 N/sample	30 N/sample	30 N/sample	30 N/sample	30 N/sample	30 N/sample	30 N/sample	630 g*	
Time	until planar	1,5 min	1,0 min	2x 0,5 min	2x 0,5 min	3,0 min	3,0 min	2 h	
Etching									
Remarks	* 1 large weight								

Legend [EN]	Legende [DE]
DPS = Diamond Plane Grinding Disc   DSS = Diamond Grinding Disc   ADAMANT = ADAMANT Diamond Grinding Disc   IDAMANT = IDAMANT Diamond Grinding and Polishing Disc	DPS = Diamant-Planschleifscheibe   DSS = Diamant-(Fein)schleifscheibe   ADAMANT = ADAMANT Diamantschleifscheibe   IDAMANT = IDAMANT-Schleifpolierscheibe
>< Countercurrent   >> Concurrent	>< Gegenlauf   >> Gleichlauf
mkd = monocrystalline Diamond, concentration 50 carat/litre pkd = polycrystalline Diamond, concentration 50 carat/litre	mkd = monokristalliner Diamant, Konzentration 50 Karat/Liter pkd = polykristalliner Diamant, Konzentration 50 Karat/Liter
Lubricant: coolTec I (Water-based)   coolTec II (Alcohol-based)   coolTec III (Oil-based)   coolTec IV (Water-based)   coolTec Orange (Anhydrous)	Schmiermittel: coolTec I (Wasserbasis)   coolTec II (Alkohobasis)   coolTec III (Ölbasis)   coolTec IV (Wasserbasis)   coolTec Orange (Wasserfrei)
Spüli: Dishwasher detergent with H <sub>2</sub> O for cleaning the samples	Spüli: Spüli mit Wasser zum Reinigen der Proben
<a href="#">You can find metallography supplies in our store (link)</a>	<a href="#">Verbrauchsmaterialien finden Sie in unserem Shop (Link)</a>